

# Lawrence Livermore National Laboratory

## Characterization of Uranium Oxyfluoride Particles for Nuclear Safeguards

6 December 2011



Ruth Kips, Ph.D.  
Staff Scientist

Lawrence Livermore National Laboratory, P. O. Box 808, Livermore, CA 94551


This work performed under the auspices of the U.S. Department of Energy by  
Lawrence Livermore National Laboratory under Contract DE-AC52-07NA27344

UCRL-PRES-496492

# The next 10 min I will talk about...

- ☐ What is **particle analysis** for **environmental sampling (ES)**?
- ☐ What are the **current trends** in nuclear safeguards?
- ☐ The uranium oxyfluoride particle project at LLNL
  - Scope of the work
  - Analytical techniques used
  - Our latest findings
  - Conclusions

# Collecting nuclear fingerprints through Environm. Sampling



Whenever nuclear material is processed, small amounts are released to the environment


Material released = characteristic of process material

Environmental sampling uses pieces of cotton to collect this material (dust, particles) **SWIPES**

Swipes analyzed in **bulk** and on a **particle basis** at IAEA's Network of Analytical Labs (NWAL)

Uranium isotopic composition is measured

Detection of undeclared nuclear materials for States w/ AP



Swipe sample

Sent to IAEA Safeguards Analytical Lab



Swipe analyzed as a whole and on a particle-by-particle basis

# The next generation of safeguards inspectors

From accountants...



...who verify declarations




...to detectives



...looking at the bigger picture

1  $\mu$ m single particle



crystallinity

surface structure

**isotopic composition**

elemental composition

shape & size

molecular information

**This information could help determine the particle's source and exposure history**

# As always, there are some challenges...


Advantage analyzing individual particles vs. swipe as a whole

□ no information lost through averaging

Natural U particles are everywhere  
**How do we find the non-natural/HEU ones?**



The particles are so incredibly small  
**How do we analyze them?**



From: Dartmouth Electron Microscope Facility/Dartmouth College

300x smaller than the eye of this ant!

Particle analysis requires **very sensitive equipment with ultra-high spatial resolution!**

# Uranium oxyfluoride particles from enrichment activities




## Uranium hexafluoride ( $\text{UF}_6$ ) for the enrichment of uranium

- Used in large quantities
- Very reactive
- Small amounts inadvertently released to the atmosphere



$\text{UF}_6$  reacts with atmospheric moisture

$\text{UO}_2\text{F}_2$  particles



## Hypothesis

1. *If we can measure the fluorine in these particles...*

- ☐ Indicator of enrichment activities (NPT)

2.  *$\text{UO}_2\text{F}_2$  is known to be affected by the environment*

- ☐ F-to-U ratio may be indicator of exposure history of collected particle





Particle analysis using NanoSIMS




# The Cameca NanoSIMS 50 at LLNL


## Secondary ion mass spectrometry at ultrahigh spatial resolution



- Surface technique
- Surface sputtered with high energy ion beam
- NanoSIMS spatial resolution down to the nanometer scale
- Ion microscope
- Produces elemental maps
- LLNL NanoSIMS only NanoSIMS used for nuclear applications




$^{55}\text{Mn}$  image of "Ada" Stardust particle - estimated 2.56 billion years old – Courtesy J. Matzel



NanoSIMS image of bacteria and flagella *Hoplonympha*  
Courtesy K. Carpenter

# The Cameca NanoSIMS 50 at LLNL



## Secondary ion mass spectrometry at ultrahigh spatial resolution





# Preparation of $\text{UO}_2\text{F}_2$ particles at the EU's Joint Research Centre, IRMM in Belgium

## Aerosol deposition chamber





## Storage in different environmental conditions at LLNL


	ETS 1	ETS 2	ThF1	ThF2
Temp	25 °C	25 °C	25 °C	25 °C
RH	<15 %	30-40 %	50 %	70-77 %
Light	diffuse	UV-light	diffuse	diffuse



# Correlation F-to-U ratio – exposure conditions



8 x 8 micron ion images of F and U



F-to-U ratio for different exposure conditions

- NanoSIMS depth profiles collected for different exposure conditions
- F-to-U ratio calculated based on difference in ion intensity
- High humidity (70 % RH) and high T (350 C) caused biggest reduction in F-to-U ratio
- F still detected after 5 years at ambient conditions

# To conclude...

## □ Particle analysis for environmental sampling

= Powerful tool for nuclear safeguards inspections

- Focus on uranium isotope analysis
- Demonstrated advantages of using complementary particle analysis techniques

## □ NanoSIMS analysis of uranium oxyfluoride particles

- F and U distribution of individual particles
- Variations in F-to-U ratio for different exposure conditions
- High humidity and high temperature accelerate F loss
- F still detected after 5 years for dry/intermediate RH conditions

**Thank you!**

# Acknowledgements

Thanks to J. Zaug and D. Grant for help with the Raman measurements and to the NNSA Office of Nonproliferation and International Security (NA-243) for their generous support

*The uranium oxyfluoride particles were prepared at IRMM expressly for the purpose of this study.*

*We also have no evidence that fugitive emissions from a uranium enrichment facility would, or would not, appear similar to the analyzed samples.*

**Thank you!**