

Lawrence Livermore National Laboratory

Multi-actinide resonance ionization mass spectrometry (RIMS) for Nuclear Forensics Applications and Rapid Response



Brett H. Isselhardt, Ph.D.

Kim B. Knight, Ph.D.

Lawrence Livermore National Laboratory, P. O. Box 808, Livermore, CA 94551

This work performed under the auspices of the U.S. Department of Energy by
Lawrence Livermore National Laboratory under Contract DE-AC52-07NA27344


LLNL-PRES-516915

Nuclear forensics: What do we mean?



- analysis of interdicted illicit nuclear or radioactive material and associated materials to provide evidence for nuclear attribution
- the **goal** is identification of forensic indicators in interdicted nuclear and radiological samples or the surrounding environment, **linking people to places, materials, events**

Measuring nuclear materials: One size does not fit all




Analytical considerations:

- Sample size?
- Spatial resolution?
- Homogenous or heterogeneous samples?
- Analytical location?
- Desired precision/accuracy?
- Time frame?



There is no one path to an answer: Multiple signatures build attribution confidence




Primary tools for characterization of nuclear materials:


- **Radiation detection / counting**
 - Alpha (α), gamma (γ), beta (β^-)
 - Selective isotopic information
- **Excitation/stimulation**
 - X-ray generation
 - Atomic absorption and emission
 - Limited analysis of complex samples and trace concentrations
- **Ion formation (mass spectrometry)**
 - Any element and isotope
 - range of concentrations, spatial resolutions
 - **Generally higher precision, but time consuming sample preparation**

Moving capabilities forward: Nuclear Forensic Research

- Improve answers?
 - Precision and Accuracy
- Faster analyses?
- Smaller sample sizes?
- Fewer people needed?
- Reliable data processing?



Reference text authored at LLNL, 2005



IAEA document by Kristo (LLNL), Smith (LLNL), Niemeyer (LLNL), and Dudder (PNNL), 2006

U and Pu-isotopes contain critical clues about material origins, processes and intent

Actinide isotope measurements demand speed, precision and accuracy, rigorous conclusions



Resonance Ionization: A Faster Path to a Better Answer

**All Mass Spectrometers
Work The Same Way:**


Make ions from the element(s) of interest
Separate ions in space or time
Detect the ions as separated masses

**Resonance ionization uses intermediate electronic states
specific to each element to form ions!**



Lasers



provide the energy
discrete photon energy
sufficient photon fluence
precision, timing, control



Almost any element can be ionized through this process!

RIMS performs chemical separation in the instrument

Resonance ionization mass spectrometry measures actinide isotopes using selective laser ionization for chemical separation during analysis




Laser Performance


Real lasers vary in:

- Wavelength
- Position
- Power
- Time


Time Variation from Pulse to Pulse:




Old laser: 14 ns FWHM




New laser: 4 ns FWHM



Wavelength Variation





Position




RIMS demonstrates discrimination between U and Pu from the same sample

- ^{238}U and ^{238}Pu would normally have to be separated chemically, prior to mass spectrometry
- **RIMS does this separation in the instrument**
- Tuned lasers ionize the desired element for analysis from a multi-actinide sample




Secondary electron image



Actinide isotope ratios have been measured accurately in metal, oxide and glass samples

Standard-bracketed sample (U970) analyses demonstrate presently achievable accuracy



CRM970 (U_3O_8)

$$^{235}\text{U}/^{238}\text{U}_{\text{cert}} = 186.77$$


$$^{235}\text{U}/^{238}\text{U}_{\text{meas}} = 185.16$$

$$(1\sigma) \pm 1.76$$


$$\chi^2 = 1.10$$

- We have measured isotope ratios over seven orders of magnitude
- Completed complementary U and Pu analysis in ~2 hours or less


Preliminary results demonstrate detection of actinides in fallout debris samples with part-per-million concentrations of actinides



Model: Predicting Ionization Behavior of Uranium



Variation of Measured Isotope Ratio as a Function of Laser Wavelength in the First Transition



Sustaining and advancing nuclear forensic capabilities

Nuclear forensics is an emerging science

By necessity, it draws on a team of people to research and innovate, provide process knowledge and materials experience, and perform technical analyses



Over 5 years, this application of the RIMS methodology has advanced from concept to demonstrated multi-element technique for measuring isotope ratios in actinides without chemical separation.

